

Print this page


The LIFE® Emergency Oxygen Units must be filled with medical grade oxygen (USP) if the indicator on the gauge is less than full, so that an adequate supply is on hand for emergency.

THE OXYGEN CYLINDER MAY BE REFILLED WITH PURE MEDICAL OXYGEN BY ANY MEDICAL/INDUSTRIAL GAS DISTRIBUTOR who uses the Compressed Gas Association (CGA) #870 PIN-INDEXED UNIVERSAL COUPLING FOR OXYGEN.

- | | |
|---|--|
| <ul style="list-style-type: none"> ▪ AGA Gases ▪ BOC Gases ▪ Airgas ▪ Praxair | <ul style="list-style-type: none"> ▪ Air Liquide Co. ▪ TriGas ▪ Air Products Co |
|---|--|

Check your Yellow Pages under WELDING EQUIPMENT GASES

PIN-INDEXED YOKE CONNECTIONS FOR MEDICAL GASES (CHROME PLATED BRASS)



Medical gas cylinders, which are used in respiration therapy, inhalation or resuscitation applications are pin-indexed to prevent interchangeability with dangerous or non compatible gases, in conformance with with CGA Standards.

